

# Advisory Notice

**Date Released: 2019-January-29**

**Posted By: Quality Control**

**Product Affected: Troponin I Improvement**

**Brief Description of Advisory: More sensitive, lower cutoff, new IFU**

To whom it may concern,

Based on market feedback and customer request, Monobind has improved its **Troponin I** (cTnI) assay to be more competitive with other systems on the market, including better low end-sensitivity and lower cut-off values. Changes are summarized below:

	<b>Previous cTnI System</b>	<b>Improved cTnI System</b>
Calibrator Value	0, 1, 3, 6, 15, 30 ng/ml	0, 0.4, 1.25, 2.5, 7.5, 20 ng/ml
Plate	Streptavidin Coated	Antibody Coated
Cut off Range/ Expected Value	< 1.3 ng/ml	< 0.4 ng/ml

All other aspects of the product remain unchanged, including part number, general procedure and total assay time.

The product page has been updated on our website and a new product Instructions for Use (**IFU**) released and can be downloaded at <http://www.monobind.com/Support-Areas/product-inserts>

The **Design Output** is available by request.

**Product Price:** There has been no price increase.

**Part Number:** There is no change to product part number

As an official of the company, I verify this information to be true and correct.

Sincerely,



Anthony Shatola  
Quality Control Director